



Centre for Cellular and Molecular Biology (CCMB) - India  
 CCMB Flow Cytometry Facility (C-406)



Design a panel and find antibodies today at  
<https://app.fluorofinder.com/ccmb>

## MoFlo XDP Configuration

LASER	LASER (NM)	BP	LP	OPTIMAL COLORS	Additional optimal fluorochrome options available on FluoroFinder.com.				
LASER	351	450 / 65	--	<b>DAPI</b> 	<b>Hoechst (33342)</b> DyLight 350 				
		670 / 30	--	<b>BUV661</b> 	<b>Hoechst (33342)</b> 				
LASER	488	529 / 28	--	<b>Alexa Fluor 488</b> 	<b>FITC</b> 	<b>GFP</b> 	ATTO 488 	BB515 	
				Cy2 	DyLight 488 	Fluorescein 	NL 493 	Vio488 	
		575 / 26	--	<b>dsRed</b> 	<b>PE</b> 				
		628 / 32	--	<b>PE-Texas Red</b> 	<b>PI</b> 	ECD 	PE-CF594 	PE-Dazzle 594 	
PE-eFluor 610 	PE-Vio615 								
670 / 30	--	<b>PE-Cy5</b> 	<b>7-AAD</b> 	PerCP-Vio700 					
785 / 62	--	<b>PE-Cy7</b> 	PE-Vio770 						
LASER	638	670 / 30	--	<b>Alexa Fluor 647</b> 	<b>APC</b> 	Alexa Fluor 633 	Alexa Fluor 635 	ATTO 655 	
				Cy5 	DyLight 633 	DyLight 650 	NL 637 		
		720 / 13	--	<b>Alexa Fluor 700</b> 	APC-R700 				
785 / 62	--	<b>APC-Cy7</b> 	APC-Alexa 750 	APC-eFluor 780 	APC-Fire 750 	APC-H7 			
				APC-Vio770 					

**Highlighted** colors are suggest by the instrument owner or manufacturer.

**Tip:** FluoroFinder.com enables interactive product selection from hundreds of antibody suppliers with machine specific spectra-viewers to guide you to an optimal panel design.

*"Reasonably-priced benchtop  
 cytometer that lives large"*



**NovoCyte**  
 Flow Cytometer

[Learn more](#)

